

OKmeter Core Blood Glucose Test Strips

Instructions for Use

Intended Use

IVD OKmeter Core Blood Glucose Test Strips are used with the OKmeter Core Blood Glucose Meters. They are intended for self-testing by people with diabetes at home and health care professionals in a clinical setting to monitor glucose concentrations in fresh capillary whole blood (drawn from the fingertips, forearm, upper arm, palm, calf or thigh), or venous whole blood. They are for testing outside the body (*in vitro* diagnostic use only). Do not use them for diagnosis of diabetes or testing on neonates.

Test Principle

Glucose in the blood samples mixes with a special chemical in the test strip and produces a small electric current. The amount of current produced changes with the amount of glucose in the blood. The glucose meter measures the strength of the current and displays the results as a blood glucose level.

Characteristics

Each test strip is plasma-calibrated[†], requiring a sample volume of 0.5 µL and taking just 5 seconds to return a test result. The test range is 20 to 600 mg/dL (or 1.1 to 33.3 mmol/L) with resolution at 1 mg/dL (0.1 mmol/L). [†]plasma-calibrated means the following :

- The traceable calibrator used is the YSI 2747 Glucose Standard, which is a NIST traceable glucose standard.
- The reference instrument used is the YSI 2300 Glucose Analyzer, which is calibrated by YSI 2747 Glucose Standard.
- The margin of error of the Calibrated YSI Glucose Analyzer measurement is 0.289 mg/dL when the blood glucose concentration is higher than 100 mg/dL, or 0.029 mg/dL when the blood glucose concentration is lower than 100 mg/dL. (1 mmol/L= 18 mg/dL).

[†]Test results produced on capillary whole blood samples by OKmeter Core System are compared with the results of the corresponding plasma samples tested by the calibrated YSI 2300 Glucose Analyzer.

Limitations of the System

OKmeter Core Blood Glucose Test Strips provide accurate results when the following considerations are observed:

- Use fresh capillary whole blood only. Do not use serum or plasma.
- The device can be used on neonates, but needs to be preset neonate test on the meter before testing. It should be tested by doctor or professionals only.
- The test strips are for single use only. Do not reuse.
- Hematocrit levels below 20% or above 60% can cause inaccurate results. Please consult your doctor if you do not know your hematocrit level.
- Dehydration may cause lower test results. If you are severely dehydrated, contact your physician immediately.
- Testing at altitudes up to 3,402 meters (or 11,161 ft) does not significantly affect results.

Storage and Handling

Please take the following precautions to ensure your OKmeter Core Blood Glucose Test Strips are effective.

- Prior to first use, ensure that the package is undamaged and closed.
- Keep the test strip vial away from sunlight and in a cool, dry place between 4~40°C (39~104°F). Do not freeze it.
- Store test strips in their original pack only. Do not put the test strips in any other container.
- Close the vial cap tightly immediately after removing a test strip. This keeps the remaining test strips fully functional right up to the expiration date.
- Use test strips immediately after removing from the package.

- Do not use test strips after the expiration date.
- The strip vial contains desiccant and indicator beads at the bottom. If the indicator beads change color from blue to pink, the test strips have been exposed to moisture and should not be used.
- To ensure accurate results, do not use expired test strips or strips with indicator beads that have changed color.
- Avoid getting dirt, food or water on the test strip. Do not handle test strips with wet hands. All parts of the test strip should be touched only with dry and clean fingers.
- Do not perform blood glucose tests at a temperature below +10°C (50°F) or above +40°C (104°F), and relative humidity between 10% and 85%.
- Make a notation of the date on the vial label when you first open it. Discard remaining test strips 180 days after first opening the vial.

Warning

- Discard used test strips and lancets responsibly according to your local regulations.
- Keep test strips vial away from children. A child could choke on the test strips. The test strips and their vial contain agents that may be harmful if swallowed. If they are swallowed, promptly see a doctor for immediate assistance.

⚠ Do not change medication based on the test results without the advice of a physician or healthcare professional.

Additional Information for Healthcare Professionals

1. Follow the infection control procedures appropriate for your facility.
2. OKmeter Core Blood Glucose Test Strips are not validated for and should not be used for testing venous blood specimens.
3. Cholesterol concentrations up to 500 mg/dL (12.9 mmol/L) or triglycerides up to 2000 mg/dL (22.6 mmol/L) do not significantly affect test results. However, glucose values in specimens beyond these levels should be interpreted with

caution.

4. Inaccurate results may occur on severely hypotensive individuals or patients in shock. Inaccurate results may also occur when individuals are in hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with glucose meters.

Interference: Reducing substances occurring in the blood naturally (uric acid, bilirubin) or from therapeutic treatments (ascorbic acid, acetaminophen) will not significantly affect results. The limiting concentration of several compounds are listed in below chart:

Compounds	Concentrations higher than the following values may cause inaccurate results
Acetaminophen	8.0 mg/dL (0.53 mmol/L)
Ascorbic Acid	5.0 mg/dL (0.28 mmol/L)
Aspirin	60 mg/dL (3.33 mmol/L)
Bilirubin	90 mg/dL (1.54 mmol/L)
Cholesterol	500 mg/dL (12.9 mmol/L)
Creatinine	5.0 mg/dL (0.44 mmol/L)
Dopamine	2.0 mg/dL (0.11 mmol/L)
EDTA	360 mg/dL (12.3 mmol/L)
Galactose	900 mg/dL (50 mmol/L)
Gentisic Acid	5.0 mg/dL (0.32 mmol/L)
Glutathione	53 mg/dL (1.72 mmol/L)
Haemoglobin	500 mg/dL (0.08 mmol/L)
Heparin	8,000 U/dL
Hydroxyurea	3.0 mg/dL (0.39 mmol/L)
Ibuprofen	50 mg/dL (2.42 mmol/L)
Icodextrin	13 mg/dL (0.01 mmol/L)
L-dopa	10 mg/dL (0.51 mmol/L)
Maltose	900 mg/dL (26.3 mmol/L)
Methyldopa	3.0 mg/dL (0.13 mmol/L)
Pralidoxime Iodide	25 mg/dL (0.94 mmol/L)
Salicylate	60 mg/dL (4.34 mmol/L)

Compounds	Concentrations higher than the following values may cause inaccurate results
Tolazamide	100 mg/dL (3.21 mmol/L)
Tolbutamide	400 mg/dL (14.8 mmol/L)
Triglycerides	2,000 mg/dL (22.6 mmol/L)
Uric Acid	8.0 mg/dL (0.48 mmol/L)
Xylose	100 mg/dL (6.66 mmol/L)

Blood Glucose Testing Procedure

See your OKmeter Core User Guide and accompanying insert for detailed illustrations for all test procedures.

Test Results

Test results are shown in milligrams of glucose per deciliter of blood (mg/dL) or in millimoles of glucose per liter (mmol/L).

The meter is capable of displaying test results in the range of 20 to 600 mg/dL (or 1.1 to 33.3 mmol/L).

Glucose levels below 60 mg/dL (3.3 mmol/L)*¹ or above 240 mg/dL (13.3 mmol/L)*² may indicate a potentially serious medical condition. If your test result is below 60 mg/dL (or 3.3 mmol/L) or above 240 mg/dL (or 13.3 mmol/L), please consult your healthcare professional immediately.

Reference values: The normal adult fasting blood glucose range for a nondiabetic person is Less than 100 mg/dL and less than 140 mg/dL up to 2 hours after meals*³.

These are expected values for people without diabetes. Users are to work with their healthcare professional to determine their target blood glucose values.

*1: Kahn, R. and Weir, G.: *Joslin's Diabetes Mellitus, 13th ed Philadelphia: Lea and Febiger (1994), 489.*

*2: Krall, L.P. and Beaser, R. S.: *Joslin Diabetes Manual. Philadelphia: Lea and Febiger (1989), 261-263.*

*3: American Diabetes Association Website (<http://www.diabetes.org>)

Inconsistent Results

If you are getting test results which are inconsistent with your state of wellness or how you feel, please do the following:

- Make sure the blood sample applied completely fills the test strip channel.
- Check that the test strips have not expired.
- Verify the performance of the meter and the test strips using the control solution.

⚠ Consult your doctor if you continue getting the same high or low results.

Quality Control (QC) Testing

Run a control test anytime you want to check the performance of the meter, the test strip or your testing technique. Only use OKmeter Control Solution. These control solution is designed specifically for use with the system. The control results should fall within the control ranges printed on the test strip bottle.

⚠ Important: the control solution range may vary with each new box of test strips. Always use the control range on the label of your current vial of test strips.

Test Strip Chemical Components

Each glucose test strip contains

- Glucose Oxidase (*Aspergillus niger*) 20 IU
- Potassium Ferricyanide 0.12 mg
- Non-reactive ingredients 1.8 mg

Performance Characteristics

Precision

Standard deviation (SD) for each glucose concentration < 100 mg/dL (5.55 mmol/L) and coefficient of variation (CV) for each glucose concentration ≥ 100 mg/dL (5.55 mmol/L) is ≤ 5.0 mg/dL (0.278 mmol/L) and ≤ 5.0%, respectively.

Intermediate precision

Control Solution Level (mg/dL)		Low (30~50)	Normal (96~144)	High (280~420)
Pooled	Mean (mg/dL)	45.5	119.3	351
	SD	1.8	3.1	10.1
	CV (%)	-	2.6%	2.9%

Repeatability

Blood Glucose (mg/dL)	30~50	51~110	111~150	151~250	251~400	
Pooled	Mean (mg/dL)	43.6	79.4	131	196	326
	SD	2.6	2.7	3.3	4.9	9.5
	CV (%)	5.9%	3.4%	2.5%	2.5%	2.9%

System Accuracy :

For glucose concentration < 100 mg/dL (5.55 mmol/L)

Within±5 mg/dL (Within±0.28mmol/L)	Within±10 mg/dL (Within±0.56mmol/L)	Within±15 mg/dL (Within±0.83mmol/L)
148/228(64.9%)	210/228(92.1%)	224/228(98.2%)

For glucose concentration ≥ 100 mg/dL (5.55 mmol/L)

Within±5%	Within±10%	Within±15%
276/444(62.2%)	393/444(88.5%)	433/444(97.5%)














For glucose concentrations between 43.4 mg/dL (2.41 mmol/L) and 446 mg/dL (27.8 mmol/L)

Within±15 mg/dL (0.83 mmol/L) or ±15%
657/672 (97.77%)

The OKmeter Core Blood Glucose Monitoring System meets the requirements for System Accuracy as stated in EN ISO 15197:2015.

User Performance

A study evaluating glucose values from fingertip, palm, forearm, upper arm, calf, and thigh capillary blood samples obtained by 100 lay persons showed the following results: 100% within ± 15 mg/dL (± 0.83 mmol/L) of the medical laboratory values at glucose concentrations below 100 mg/dL (5.55 mmol/L), and fingertip 98.4% / palm 96.7% / forearm 100% / upperarm 98.4% / calf 95.1% / thigh 96.7% within ±15% of the medical laboratory values at glucose concentrations at or above 100mg/dL (5.55 mmol/L).

Symbol	Description
	For <i>in vitro</i> diagnostic use.
	Please consult instructions for use
	Do not re-use
	Lot Number
	Caution, consult accompanying document
	Keep dry
	Keep away from sunlight
	Humidity limitation
	Temperature limitation / Store at
	Use by / Expiry date
	Manufacturer
	EU representative
	This product fulfils the requirements of Directive 98/79/EC on <i>in vitro</i> diagnostic medical device.

 OK Biotech Co., Ltd. 

No. 91, Sec. 2, Gongdao 5th Road, 30070 Hsinchu City, TAIWAN

 MDSS GmbH

Schiffgraben 41, 30175 Hannover, GERMANY

Made in Taiwan